

J2 MEDICAL LAUNCHES HIGHLY-ANTICIPATED AKUCAL EVO™ Jarrid Henderson declares a design EVOLUTION, “100% custom functionality”

PITTSBURGH, July 2009 - Since its inception, the Akucal™ has been the preferred calibration unit for all orthopaedic preoperative planning in the world's largest and most progressive hospitals, clinics, and private practice facilities. The brand-new Akucal Evolution™ is the result of tremendous team-work and communication between the J2 development team and its valued clients around the world. Jarrid Henderson, co-owner of J2 Medical reflected, “This is the result of problem solving becoming a personal vendetta.”

Efficiency: Time is invaluable; luckily, the EVO™ saves time for techs. Whether a tech needs it suctioned to a table, wall bucky or textured surface, they will be able to get the EVO™ in the desired position within seconds. The quick release suction base option allows the tech to easily release the EVO™ and move on to take the next image. The simplistic nature of the EVO™ allows it to have a full line of replaceable parts; making troubleshooting a breeze—saving valuable time and money. The indicator clip references the sphere size in every image via the solid numbers attached to the Ball Claw™ coupler. The EVO™ spheres are solid, unlike its competitors, which use hollow spheres that render a dull image. In addition, J2 Medical is releasing Smoothees™, which make any textured surface suctionable.

Durability: Henderson eliminated all plastic base parts and re-designed the EVO™ with the highest quality materials available. The EVO™ boasts machined stainless steel and anodized aluminum bases, a fully-threaded arm to base function, and a patent pending advanced Ball Claw™ coupler. The new EVO™ Flat Base is solid stainless steel and offers the smallest footprint available. The EVO™ is custom-built for very high use facilities that use the Akucal™ all day, every day.

Sanitation: The EVO™ was also designed with sanitation in mind. Because the EVO's™ parts are interchangeable and easily removed, it can be fully-cleaned and autoclaved. The entire EVO™ can be sanitized from top to bottom without wasting precious time. “This was a very important issue during the design phase—having the ability to thoroughly clean the Akucal in minimal time was one of the most vital requests of our customers,” Henderson commented.

J2 Medical Solutions, established in 2007 and located in Pittsburgh, PA is a world leader in Orthopaedic devices; and a leading developer of radiology accessories.

To learn how you can save \$100 on the purchase of your EVO™ and other J2 products please call 412- 573-6116 or visit www.j2medical.com .